

Product datasheet for **GA116114**

Human Matriptase 2 (TMPRSS6) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMPRSS6 CRISPRa kit - CRISPR gene activation of human transmembrane serine protease 6
Symbol:	TMPRSS6
UniProt ID:	Q8IU80
Synonyms:	IRIDA
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)