

## Product datasheet for **GA116118**

### Human WDR69 (DAW1) activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                            |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                               |
| Description:  | DAW1 CRISPRa kit - CRISPR gene activation of human dynein assembly factor with WD repeats 1                                                                                                                                                                                                                |
| Symbol:       | DAW1                                                                                                                                                                                                                                                                                                       |
| Components:   | GA116118G1, WDR69 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA116118G2, WDR69 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA116118G3, WDR69 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q8N136</a>                                                                                                                                                                                                                                                                                     |
| Synonyms:     | ODA16; WDR69                                                                                                                                                                                                                                                                                               |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                    |

#### Product images:



[View online »](#)