

Product datasheet for **GA116121**

Human GPR113 (ADGRF3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ADGRF3 CRISPRa kit - CRISPR gene activation of human adhesion G protein-coupled receptor F3
Symbol:	ADGRF3
UniProt ID:	Q8IZF5
Synonyms:	GPR113; PGR23
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)