


Product datasheet for **GA116124**

Human **TOGARAM2** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	TOGARAM2 CRISPRa kit - CRISPR gene activation of human TOG array regulator of axonemal microtubules 2
Symbol:	TOGARAM2
Components:	GA116124G1, TOGARAM2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116124G2, TOGARAM2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116124G3, TOGARAM2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6ZUX3</u>
Synonyms:	family with sequence similarity 179, member A; RGD1310677
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)