


Product datasheet for **GA116130**

Human PARP15 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PARP15 CRISPRa kit - CRISPR gene activation of human poly(ADP-ribose) polymerase family member 15
Symbol:	PARP15
Components:	GA116130G1, PARP15 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116130G2, PARP15 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116130G3, PARP15 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q460N3
Synonyms:	ARTD7; BAL3; pART7
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)