

Product datasheet for **GA116147**

Human TIGD2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TIGD2 CRISPRa kit - CRISPR gene activation of human tiger transposable element derived 2
Symbol:	TIGD2
Components:	GA116147G1, TIGD2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116147G2, TIGD2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116147G3, TIGD2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q4W5G0</u>
Synonyms:	HEL106
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)