

Product datasheet for **GA116158**

Human ZNF366 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF366 CRISPRa kit - CRISPR gene activation of human zinc finger protein 366
Symbol:	ZNF366
Components:	GA116158G1, ZNF366 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116158G2, ZNF366 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116158G3, ZNF366 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8N895
Synonyms:	DCSCRIPT
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)