

Product datasheet for **GA116160**

Human PRSS35 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PRSS35 CRISPRa kit - CRISPR gene activation of human serine protease 35
Symbol:	PRSS35
UniProt ID:	<u>Q8N3Z0</u>
Synonyms:	C6orf158; dj223E3.1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)