


Product datasheet for **GA116165**

Human C6orf118 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C6orf118 CRISPRa kit - CRISPR gene activation of human chromosome 6 open reading frame 118
Symbol:	C6orf118
Components:	GA116165G1, C6orf118 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116165G2, C6orf118 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116165G3, C6orf118 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5T5N4</u>
Synonyms:	bA85G2.1; dj416F21.2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)