

Product datasheet for **GA116174**

Human PKD1L1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PKD1L1 CRISPRa kit - CRISPR gene activation of human polycystin 1 like 1, transient receptor potential channel interacting
Symbol:	PKD1L1
Components:	GA116174G1, PKD1L1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116174G2, PKD1L1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116174G3, PKD1L1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8TDX9</u>
Synonyms:	PRO19563
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)