

Product datasheet for **GA116183**

Human **SNX31** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SNX31 CRISPRa kit - CRISPR gene activation of human sorting nexin 31
Symbol:	SNX31
Components:	GA116183G1, SNX31 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116183G2, SNX31 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116183G3, SNX31 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8N9S9</u>
Synonyms:	MGC39715
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)