

Product datasheet for **GA116184**

Human TMEM64 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM64 CRISPRa kit - CRISPR gene activation of human transmembrane protein 64
Symbol:	TMEM64
Components:	GA116184G1, TMEM64 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116184G2, TMEM64 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116184G3, TMEM64 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6YI46</u>
Synonyms:	DKFZp762C1112
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)