

Product datasheet for **GA116188**

Human **KCNV2** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	KCNV2 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel modifier subfamily V member 2
Symbol:	KCNV2
UniProt ID:	Q8TDN2
Synonyms:	Kv8.2; KV11.1; RCD3B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)