


Product datasheet for **GA116193**

Human ZNF883 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF883 CRISPRa kit - CRISPR gene activation of human zinc finger protein 883
Symbol:	ZNF883
Components:	GA116193G1, ZNF883 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116193G2, ZNF883 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116193G3, ZNF883 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P0CG24</u>
Synonyms:	RP11-349F12.1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)