

Product datasheet for **GA116195**

Human TENT5D activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TENT5D CRISPRa kit - CRISPR gene activation of human terminal nucleotidyltransferase 5D
Symbol:	TENT5D
Components:	GA116195G1, TENT5D gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116195G2, TENT5D gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116195G3, TENT5D gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8NEK8</u>
Synonyms:	CT1.26; CT112
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)