

Product datasheet for GA116201

Human ARX activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | ARX CRISPRa kit - CRISPR gene activation of human aristaless related homeobox |
| Symbol: | ARX |
| UniProt ID: | Q96QS3 |
| Synonyms: | CT121; EIEE1; ISSX; MRX29; MRX32; MRX33; MRX36; MRX38; MRX43; MRX54; MRX76; MRX87; MRX51; PRTS |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)