

Product datasheet for **GA116211**

Human HTR3C activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HTR3C CRISPRa kit - CRISPR gene activation of human 5-hydroxytryptamine receptor 3C
Symbol:	HTR3C
Components:	GA116211G1, HTR3C gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116211G2, HTR3C gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116211G3, HTR3C gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8WXA8
Synonyms:	5-hydroxytryptamine (serotonin) receptor 3, family member C; 5-hydroxytryptamine receptor 3 subunit C
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)