

Product datasheet for **GA116227**

Human **KCN G3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	KCN G3 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel modifier subfamily G member 3
Symbol:	KCN G3
UniProt ID:	Q8TAE7
Synonyms:	KV6.3; KV10.1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)