


Product datasheet for **GA116228**

Human Phostensin (PPP1R18) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PPP1R18 CRISPRa kit - CRISPR gene activation of human protein phosphatase 1 regulatory subunit 18
Symbol:	PPP1R18
Components:	GA116228G1, Phostensin gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116228G2, Phostensin gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116228G3, Phostensin gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6NYC8</u>
Synonyms:	HKMT1098; KIAA1949
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)