

Product datasheet for **GA116233**

Human ZNF384 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF384 CRISPRa kit - CRISPR gene activation of human zinc finger protein 384
Symbol:	ZNF384
UniProt ID:	Q8TF68
Synonyms:	CAGH1; CAGH1A; CIZ; ERDA2; FLJ59043; NMP4; NP; TNRC1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)