

Product datasheet for **GA116250**

Human PGAM5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PGAM5 CRISPRa kit - CRISPR gene activation of human PGAM family member 5, mitochondrial serine/threonine protein phosphatase
Symbol:	PGAM5
UniProt ID:	<u>Q96HS1</u>
Synonyms:	BXLBV68
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)