

Product datasheet for **GA116264**

Human MPZL3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MPZL3 CRISPRa kit - CRISPR gene activation of human myelin protein zero like 3
Symbol:	MPZL3
Components:	GA116264G1, MPZL3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116264G2, MPZL3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116264G3, MPZL3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6UWV2</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)