

Product datasheet for **GA116269**

Human DNAH10 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DNAH10 CRISPRa kit - CRISPR gene activation of human dynein axonemal heavy chain 10
Symbol:	DNAH10
Components:	GA116269G1, DNAH10 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116269G2, DNAH10 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116269G3, DNAH10 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8IVF4</u>
Synonyms:	DNAH10 variant protein; dynein, axonemal, heavy chain 10; dynein, axonemal, heavy polypeptide 10; FLJ38262; FLJ38262, FLJ43486, FLJ43808, KIAA2017; FLJ43486; FLJ43808; KIAA2017
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

