

Product datasheet for **GA116298**

Human ZNF778 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF778 CRISPRa kit - CRISPR gene activation of human zinc finger protein 778
Symbol:	ZNF778
Components:	GA116298G1, ZNF778 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116298G2, ZNF778 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116298G3, ZNF778 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96MU6
Synonyms:	FLJ31875; MGC150573
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)