


Product datasheet for **GA116300**

Human WDR90 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WDR90 CRISPRa kit - CRISPR gene activation of human WD repeat domain 90
Symbol:	WDR90
Components:	GA116300G1, WDR90 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116300G2, WDR90 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116300G3, WDR90 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96KV7</u>
Synonyms:	C16orf15; C16orf16; C16orf17; C16orf18; C16orf19
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)