

Product datasheet for **GA116307**

Human **TTC21A** activation kit by **CRISPRa**

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | TTC21A CRISPRa kit - CRISPR gene activation of human tetratricopeptide repeat domain 21A |
| Symbol: | TTC21A |
| UniProt ID: | <u>Q8NDW8</u> |
| Synonyms: | IFT139A; STI2 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)