

## Product datasheet for **GA116316**

### Human U2AF1L3 (U2AF1L4) activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | U2AF1L4 CRISPRa kit - CRISPR gene activation of human U2 small nuclear RNA auxiliary factor 1 like 4 |
| Symbol:       | U2AF1L4  |
| UniProt ID:   | <a href="#">Q8WU68</a>   |
| Synonyms:     | U2AF1-RS3; U2AF1L3; U2AF1L3V1; U2AF1RS3; U2af26  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                              |

#### Product images:



[View online >](#)