


Product datasheet for **GA116320**

Human LCE4A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LCE4A CRISPRa kit - CRISPR gene activation of human late cornified envelope 4A
Symbol:	LCE4A
Components:	GA116320G1, LCE4A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116320G2, LCE4A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116320G3, LCE4A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5TA78</u>
Synonyms:	LEP8; SPRL4A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)