


Product datasheet for **GA116330**

Human **CC2D1B** activation kit by **CRISPRa**

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | CC2D1B CRISPRa kit - CRISPR gene activation of human coiled-coil and C2 domain containing 1B |
| Symbol: | CC2D1B |
| Components: | GA116330G1, CC2D1B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116330G2, CC2D1B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116330G3, CC2D1B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID: | <u>Q5T0F9</u> |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online »](#)