

Product datasheet for **GA116352**

Human TSC21 (TEX37) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TEX37 CRISPRa kit - CRISPR gene activation of human testis expressed 37
Symbol:	TEX37
UniProt ID:	<u>Q96LM6</u>
Synonyms:	C2orf51; TSC21
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)