

## Product datasheet for **GA116355**

### Human PIP5K3 (PIKFYVE) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PIKFYVE CRISPRa kit - CRISPR gene activation of human phosphoinositide kinase, FYVE-type zinc finger containing
Symbol:	PIKFYVE
UniProt ID:	<a href="#">Q9Y2I7</a>
Synonyms:	CFD; KIAA0981; MGC40423; PIKfyve; PIP5K
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)