

## Product datasheet for **GA116364**

### Human LIPH activation kit by CRISPRa

#### Product data:

|               |                                                                         |
|---------------|-------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                            |
| Description:  | LIPH CRISPRa kit - CRISPR gene activation of human lipase H             |
| Symbol:       | LIPH                                                                    |
| UniProt ID:   | <u>Q8WWY8</u>                                                           |
| Synonyms:     | AH; ARWH2; HYPT7; LAH2; LPDLR; mPA-PLA1; PLA1B                          |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

#### Product images:



[View online >](#)