

## Product datasheet for **GA116376**

### Human **FAM18B2 (TVP23C)** activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TVP23C CRISPRa kit - CRISPR gene activation of human trans-golgi network vesicle protein 23 homolog C
Symbol:	TVP23C
UniProt ID:	<u><a href="#">Q96ET8</a></u>
Synonyms:	FAM18B2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)