

## Product datasheet for **GA116380**

### Human SH3D20 (ARHGAP27) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	ARHGAP27 CRISPRa kit - CRISPR gene activation of human Rho GTPase activating protein 27
Symbol:	ARHGAP27
UniProt ID:	<u><a href="#">Q6ZUM4</a></u>
Synonyms:	CAMGAP1; FLJ35572; FLJ43547; MGC120624
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)