

Product datasheet for **GA116381**

Human ZNF385C activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF385C CRISPRa kit - CRISPR gene activation of human zinc finger protein 385C
Symbol:	ZNF385C
UniProt ID:	Q66K41
Synonyms:	zinc finger protein 385C
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)