

Product datasheet for **GA116383**

Human C17orf108 (LYRM9) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LYRM9 CRISPRa kit - CRISPR gene activation of human LYR motif containing 9
Symbol:	LYRM9
Components:	GA116383G1, C17orf108 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116383G2, C17orf108 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116383G3, C17orf108 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A8MSI8</u>
Synonyms:	C17orf108; HSD24
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)