


Product datasheet for **GA116387**

Human **LRRC45** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	LRRC45 CRISPRa kit - CRISPR gene activation of human leucine rich repeat containing 45
Symbol:	LRRC45
Components:	GA116387G1, LRRC45 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116387G2, LRRC45 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116387G3, LRRC45 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96CN5</u>
Synonyms:	MGC20806
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)