

Product datasheet for **GA116388**

Human SLC39A11 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLC39A11 CRISPRa kit - CRISPR gene activation of human solute carrier family 39 member 11
Symbol:	SLC39A11
Components:	GA116388G1, SLC39A11 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116388G2, SLC39A11 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116388G3, SLC39A11 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8N1S5</u>
Synonyms:	C17orf26; ZIP11
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)