

Product datasheet for **GA116392**

Human SPNS3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SPNS3 CRISPRa kit - CRISPR gene activation of human sphingolipid transporter 3 (putative)
Symbol:	SPNS3
Components:	GA116392G1, SPNS3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116392G2, SPNS3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116392G3, SPNS3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6ZMD2</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

