

## Product datasheet for **GA116403**

### Human TIGIT activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TIGIT CRISPRa kit - CRISPR gene activation of human T cell immunoreceptor with Ig and ITIM domains
Symbol:	TIGIT
UniProt ID:	<a href="#">Q495A1</a>
Synonyms:	VSIG9; VSTM3; WUCAM
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)