

Product datasheet for **GA116406**

Human TIGD4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TIGD4 CRISPRa kit - CRISPR gene activation of human tiger transposable element derived 4
Symbol:	TIGD4
Components:	GA116406G1, TIGD4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116406G2, TIGD4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116406G3, TIGD4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8IY51
Synonyms:	MGC43837
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)