

Product datasheet for **GA116407**

Human TMEM154 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM154 CRISPRa kit - CRISPR gene activation of human transmembrane protein 154
Symbol:	TMEM154
Components:	GA116407G1, TMEM154 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116407G2, TMEM154 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116407G3, TMEM154 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6P9G4</u>
Synonyms:	FLJ32028
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

