

Product datasheet for **GA116408**

Human C4orf34 (SMIM14) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SMIM14 CRISPRa kit - CRISPR gene activation of human small integral membrane protein 14
Symbol:	SMIM14
Components:	GA116408G1, C4orf34 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116408G2, C4orf34 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116408G3, C4orf34 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96QK8</u>
Synonyms:	C4orf34
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)