

Product datasheet for **GA116416**

Human RANBP3L activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	RANBP3L CRISPRa kit - CRISPR gene activation of human RAN binding protein 3 like
Symbol:	RANBP3L
Components:	GA116416G1, RANBP3L gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116416G2, RANBP3L gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116416G3, RANBP3L gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q86VV4</u>
Synonyms:	FLJ25422
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)