

Product datasheet for **GA116420**

Human **STK32A** activation kit by **CRISPRa**

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | STK32A CRISPRa kit - CRISPR gene activation of human serine/threonine kinase 32A |
| Symbol: | STK32A |
| UniProt ID: | Q8WU08 |
| Synonyms: | YANK1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)