

Product datasheet for **GA116435**

Human **LGI3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	LGI3 CRISPRa kit - CRISPR gene activation of human leucine rich repeat LGL family member 3
Symbol:	LGI3
Components:	GA116435G1, LGI3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116435G2, LGI3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116435G3, LGI3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8N145</u>
Synonyms:	LGIL4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)