

Product datasheet for **GA116436**

Human C9ORF91 (TMEM268) activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | TMEM268 CRISPRa kit - CRISPR gene activation of human transmembrane protein 268 |
| Symbol: | TMEM268 |
| UniProt ID: | Q5VZI3 |
| Synonyms: | RP11-402G3.2 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)