


Product datasheet for **GA116445**

Human SLC25A43 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLC25A43 CRISPRa kit - CRISPR gene activation of human solute carrier family 25 member 43
Symbol:	SLC25A43
Components:	GA116445G1, SLC25A43 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116445G2, SLC25A43 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116445G3, SLC25A43 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8WUT9</u>
Synonyms:	mitochondrial solute carrier protein; OTTHUMP00000023899; solute carrier family 25, member 43
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)