

Product datasheet for **GA116446**

Human RTL3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	RTL3 CRISPRa kit - CRISPR gene activation of human retrotransposon Gag like 3
Symbol:	RTL3
UniProt ID:	<u>Q8N8U3</u>
Synonyms:	Mar3; Mart3; ZHC5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)