

Product datasheet for **GA116455**

Human LASS3 (CERS3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CERS3 CRISPRa kit - CRISPR gene activation of human ceramide synthase 3
Symbol:	CERS3
UniProt ID:	Q8IU89
Synonyms:	ARCI9; LASS3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)